Optimum Growth™ Flasks – Brief Instructions

INFORS **HT**

Culture medium

Commercially available culture media are suitable. **In case that they do not contain surfactants, you must add Pluronic (Life Technologies, 0.05 g/L to 0.5 g/L) or Simethicone (Sigma, 1 %).** This reduces shear stress on the cells and foaming. **Glucose feed recommendation:** start with 80 % volume, add 10 % on day 2 and day 6. **Culture duration:** up to 20 days.

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Shaking speed

Optimum Growth™ Flasks can provide increased aeration over standard shake flasks, but usual shaking speeds of 80–120 min⁻¹ (25 mm shaking throw) are not enough. **Use 70 min⁻¹ for Hybridoma cells, 150 min⁻¹ for CHO, HEK and insect cells (25 mm shaking throw).** If you have used spinner flasks and roller bottles so far, your cells may require one or two passages at lower shaking speeds to adapt to the movement. INFORS HT recommends to use Optimum Growth™ Flasks only on 25 mm shaking diameter. 50 mm: reduce min⁻¹ by 25–40 %.

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Filling volume

Optimum Growth™ Flasks can be used with up to 50 % (with feeding up to 60 %) of the nominal volume. **Maximum gas transfer is achieved with about 30 % for Hybridomas, 50 % with CHO, HEK and insect cells.** Perform the first experiments with the optimum filling volume to see how much biomass you can get. Then, use the maximum filling volume, which may lead to less biomass/protein yield.

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