



Benchmark

Custom Designed
Clinical & Production
Freeze Dryers


designed to help you build your business



document # 234 rev 1

Custom Designed to Meet Your Requirements

- Pilot
 - Small Scale
- Sterile or Non Sterile
 - Bulk (to 320 sq ft)
- Stoppering (to 240 sq ft)
 - API Manufacturing



SP Scientific Benchmark Lyophilizers - Shelf Area & Clearance											
Model		1000	2000	3000	3500	4000	4500	5000	5500	6000	6500
Shelf Size	mm	410 x 610	510 x 760	610 x 915	610x 1220	610x 915	610x 1220	915 x 1220	1220x 1220	915x 1220	1220x 1220
	Inches	16 x 24	20 x 30	24 x 36	24 x 48	24 x 36	24 x 48	36 x 48	48 x 48	36 x 48	48 x 48
2	Area m ² (ft ²)	0.49 (5.3)	0.77 (8.33)								
	Clearance (mm)	163	215								
3	Area m ² (ft ²)	0.74 (8.0)	1.16 (12.0)	1.67 (18.0)							
	Clearance (mm)	108	139	167							
4	Area m ² (ft ²)	0.99 (10.47)	1.55 (16.67)	2.23 (24.0)	2.97 (32.0)	2.28 (24.0)	2.97 (32.0)	4.48 (48.0)	5.94 (64.0)		
	Clearance (mm)	78	101	122	122	225	225	270	270		
5	Area m ² (ft ²)	1.24 (13.33)	1.93 (20.8)	2.79 (30.0)	3.71 (40.0)	2.79 (30.0)	3.71 (40.0)	5.61 (60.0)	7.43 (80.0)		
	Clearance (mm)	60	78	95	95	175	175	210	210		
6	Area m ² (ft ²)	1.5 (16.0)	2.32 (25.0)	3.34 (36.0)	4.46 (48.0)	3.34 (36.0)	4.46 (48.0)	6.72 (72.0)	8.91 (96.0)	6.72 (72.0)	8.91 (96.0)
	Clearance (mm)	48	63	77	77	145	145	170	170	255	255
7	Area m ² (ft ²)		2.7 (29.17)	3.9 (42.0)	5.2 (56.0)	3.9 (42.0)	5.2 (56.0)	7.85 (84.0)	10.4 (112.0)	7.85 (84.0)	10.4 (112.0)
	Clearance (mm)		52	64	64	120	120	145	145	215	215
8	Area m ² (ft ²)			4.46 (48.0)	5.94 (64.0)	4.46 (48.0)	5.94 (64.0)	8.97 (96.0)	11.89 (128.0)	8.97 (96.0)	11.89 (128.0)
	Clearance (mm)			54	54	105	105	125	125	185	185
9	Area m ² (ft ²)					5.05 (54.0)	6.68 (72.0)	10.09 (108.0)	13.37 (144.0)	10.09 (108.0)	13.37 (144.0)
	Clearance (mm)					92	92	110	110	163	163
10	Area m ² (ft ²)					5.61 (60.0)	7.43 (80.0)	11.21 (120.0)	14.86 (160.0)	11.21 (120.0)	14.86 (160.0)
	Clearance (mm)					82	82	95	95	145	145
11	Area m ² (ft ²)					6.17 (66.0)	8.17 (88.0)	12.34 (132.0)	16.35 (176.0)	12.34 (132.0)	16.35 (176.0)
	Clearance (mm)					72	72	85	85	130	130
12	Area m ² (ft ²)					6.72 (72.0)	8.91 (96.0)	13.45 (144.0)	17.83 (192.0)	13.45 (144.0)	17.83 (192.0)
	Clearance (mm)					65	65	75	75	115	115
13	Area m ² (ft ²)					7.28 (78.0)	9.6 (104.0)	14.58 (156.0)	19.32 (208.0)	14.58 (156.0)	19.32 (208.0)
	Clearance (mm)					60	60	70	70	105	105
14	Area m ² (ft ²)							15.7 (168.0)	20.81 (224.0)	15.7 (168.0)	20.81 (224.0)
	Clearance (mm)							63	63	100	100
15	Area m ² (ft ²)									16.82 (180.0)	22.29 (240.0)
	Clearance (mm)									90	90



6000	6500
915x 1220	1220x 1220
36 x 48	48 x 48

15	Area m ² (ft ²)	16.82 (180.0)	22.29 (240.0)
	Clearance (mm)	90	90
16	Area m ² (ft ²)	17.34 (192.0)	23.78 (256.0)
	Clearance (mm)	85	85
17	Area m ² (ft ²)	19.06 (204.0)	25.26 (272.0)
	Clearance (mm)	75	75
18	Area m ² (ft ²)	20.18 (216.0)	26.75 (288.0)
	Clearance (mm)	2,872	72
19	Area m ² (ft ²)	21.3 (228.0)	28.24 (304.0)
	Clearance (mm)	65	65
20	Area m ² (ft ²)	22.42 (240.0)	29.72 (320.0)
	Clearance (mm)	63	63

Benchmark with SIP



Benchmark 5000




Benchmark 3000 with CIP



Benchmark 5000



Convenient "service side" layout
of clean room installation.
Designed for easy accessibility to
major components and maintenance points.



Benchmark 5000



alternate view of "service side"











Benchmark Design Flexibility System Installed in Isolator




Isolator Interface











Benchmark 5000 Internal Condenser

**When space saving is
paramount and chamber
isolation is not required**



VirTis
SP SCIENTIFIC

Benchmark 5000 in Factory Acceptance Testing

**Factory supplied FAT
documentation helps
shorten commission
time and validates
system performance**



Benchmarks Design Flexibility to Meet any Need

- **External & Internal Condensers.**
 - **Steam sterilizable (SIP) option**
- **Clean-in-place (CIP) options**




Shelf Array Performance

Shelf to shelf $\pm 1.0^{\circ}\text{C}$ at common measuring points under steady state conditions.

AISI type 316L stainless steel.

Ra25 standard Chamber and RA 32 condenser interior finish; an Ra 12 mirror finish is available as an options as well as electro-polishing.



Shelves


Shelf stack temperature uniformity - heat transfer fluid circuit designed in *parallel* (rather than in sequence)

Stainless braided Teflon hoses transport the fluid in and out of the shelves.


A compensating radiant shelf is included ensuring all product is under the same conditions

Full shelf compression is achieved in all Benchmark stoppering models.

If room height is a problem we can reduce height of cylinder, and the shelves will have "daylight" between them at full stroke




Bottom Up Stoppering



Bulk Fixed Shelf Stack

(watch the load height)



ControLyo™ Nucleation On-Demand Technology

- Load vials onto shelf of Lyophilizer at 4°C and equilibrate
- Pressurize system with inert gas (Nitrogen or Argon)
- Ramp temperature down to -4°C and begin equilibration for 45'
- Depressurize to induce nucleation

All Steps are Programmed Into the Software

- Continue to ramp down to -45°C to complete freezing step



Lynos Software with ControlLyTM Nucleation On-Demand Technology

FTS Systems **SP SCIENTIFIC**
Making Time for Science

SYSTEM ADMINISTRATOR 12/28/2010 11:56:25

Freeze Dry Pressurization Active

Batch ID: **Ballistic** Recipe Name: **OC Run** Product Name: **Chris B.**

Freezing/Thermal Treatment:

Step	1	2	3	4	5	6	7	8	9	10	11	12
Temperature SP	4.0	4.0	-45.0									
Ramp Time	0	16	90									
Hold Time	60	90	1440									

Control With CM Guaranteed Soak Product Driven Mode Product Temperature

Nucleation Control Pressurize After Step **1** Depressurize After Step **2** Shelf Load

Shelf Temp SP **-70.0** °C Hold Time **60** min Condenser SP **-30.0** °C Initial Vacuum SP **100** mTorr

Primary:

Step	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17
Temperature SP	70.0	60.0	50.0	25.0	40.0	20.0	10.0	0.0	40.0	20.0	20.0	30.0	40.0	20.0	40.0	-10.0	40.0
Ramp Time	0	60	0	0	300	20	60	0	0	20	120	480	20	0	0	0	120
Hold Time	60	20	30	60	30	30	20	30	60	30	60	30	60	120	60	60	60
Vacuum SP	20	800	500	500	100	400	200	300	50	0	300	100	300	125	75	150	400
Vacuum Ramp Time	0	30	30	0	0	20	0	30	30	0	0	0	60	0	0	0	120
Pressure Rise SP	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
PVG/CM DIFF	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0

Secondary:

Product Storage: Temperature SP **20.0** °C Vacuum SP **0** mTorr


Backfill/Stoppering: Backfill mBar Stoppering End Cycle

PVG/CM DIFF. Pressure Rise Control Action

Retest Time **2** min Closure Time **60** sec


Synoptic Freeze Dry SMART Leak Test Function Test Semi-Auto Edit Recipe

Typical Freezing Curve of Water



Uncontrolled Nucleation

SPS_HistoricalTrend.grf




PROD PROBE 09	-13.7
PROD PROBE 10	-12.9
PROD PROBE 11	-13.5
PROD PROBE 12	-15.9
PROD PROBE 13	-11.9
PROD PROBE 14	-14.2
PROD PROBE 15	-16.4
PROD PROBE 16	-14.4
THERMOCOUPLE PROBE AVERAGE	-13.8

Temperature/Vacuum

Controlled Nucleation


Nucleation Actuated




12:09:38 12/15/2010 12:10:58 12/15/2010 12:12:18 12/15/2010 12:13:38 12/15/2010

Time

PROD PROBE 09	-4.4
PROD PROBE 10	-4.3
PROD PROBE 11	-4.2
PROD PROBE 12	-4.5
PROD PROBE 13	-4.6
PROD PROBE 14	-4.0
PROD PROBE 15	-4.7
PROD PROBE 16	-4.9
THERMOCOUPLE PROBE AVERAGE	-4.1




Improved Uniformity




Primary Drying: Sublimation Front

"DRIED" PRODUCT
SUBLIMATION FRONT
FROZEN PRODUCT



(HEATING)






Experimental Demonstration of Benefits of Controlled Nucleation

- Increase in pore size and reduction in dry layer resistance
- Faster primary drying
- Reduced freezing stress on biologicals
- Reduced reconstitution times
- Drying time optimization
- Improved cake appearance
- Reduced protein aggregation

Sever, R. 2010. Controlled Nucleation in Lyophilization. Effects on Process and Products. CPPR Freeze Drying Pharmaceuticals and Biologicals Conference. Oct. 1, 2010. Garmisch, Germany.



New Poly-Bottom Trays (Dissassembled)




**Poly bottom captured
and upside down**




**Partially assembled with
sheet ready to be captured
between 'bands'**

**Detail showing poly
between 'bands'**




Fully assembled tray



Optional Auto Locking Doors


Eliminates potential problems caused by manual latching

All latching occurs simultaneously and saves time with door closure



VirTis
SP SCIENTIFIC

Auto locking door details



(FRONT VIEW)

- Automatic locking**
- Separation between clean and mechanical spaces**
- Hydraulic and pneumatic versions**
- SIP ASME code compatible**
- Compatible with cylindrical or rectangular chambers**



"T" Handle Doors




**Available as a cost saving for SIP units
Simplest technology**



VirTis
SP SCIENTIFIC

Auto-locking Constant Level Loading Door



Door & Isolator Interface Detail



Full Dress Panels





Dressed



VirTis
SP SCIENTIFIC


Isolation Valve & Actuator Valve Sized to Avoid Vapor Restriction



Internal vs External Condensers

Internal Condenser

- Cost efficiency
- Space saving
- Short vapor path
- Accessibility




External Condenser

- Barometric control is possible (iso valve)
- Faster turn around
- Less chance of oil back-streaming
- Temperatures more uniform
- Higher ice capacities





VirTis  SP SCIENTIFIC



Condenser Load Testing at Factory with Customer Verification

External Vertical Condenser



Top View



Refrigeration Options



SL- Single compressor direct expansion

-60 C shelves –70C Condenser

EL – Cascade system direct expansion

-70 C Shelves –80C Condenser


LL - Dual compressors in parallel direct expansion

-60 C Shelves –75C Condenser

RSL –Redundant system fluid condenser

-60C shelf –65 C Condenser

Upgrades to Screw compressors




Typical SL Synoptic (External)

Basic Components:

White – Vacuum


Yellow – Heat Transfer

Blue – Refrigeration



Carlyle Reciprocating Compressors


- Worldwide availability
- Rugged & dependable
- Fiscally responsible



Refrigeration Gauges




Heat Transfer Pump





Virtis  SP SCIENTIFIC

Leybold Vacuum Pumps





Vacuum Pumps with Roots Blowers

Steam in Place (SIP)



Injectibles SIP
Longer Lead Time
ASME Pressure vessel
Orbital Welding
\$Paperwork\$





Pressure and Vacuum sensors mounted on Triclover connections




Heated Head Capacitance Manometer


Ceramic Thermocouple Jacks




0.2 Micron Hydrophobic N2 Filter



Clean Diaphragm Valves



Clean-In-Place (CIP)



- Cleaning in place includes spray nozzles in the chamber and condenser
- Used to remove cyto-toxic products or insoluble residues



Clean in Place



06/15/2006 07:20:18

Programming CIP Phase **1**

Single Pass or Recirculated Wash **RECIRC**

Use CIP Skid? **No**

	Time	A1	A2	B1	B2	UP	DOWN
Valve Selection Step 1	1 Min.	On	Off	Off	Off	Off	Off
Valve Selection Step 2	2 Min.	Off	On	Off	Off	On	Off
Valve Selection Step 3	3 Min.	Off	Off	On	Off	On	On
Valve Selection Step 4	4 Min.	Off	Off	Off	On	Off	Off
Valve Selection Step 5	5 Min.	On	On	On	On	Off	On
Valve Selection Step 6	6 Min.	Off	Off	Off	On	Off	Off
Valve Selection Step 7	7 Min.	Off	Off	On	Off	Off	Off
Valve Selection Step 8	8 Min.	Off	On	Off	Off	Off	Off

CV Enabled if C1 Open **NO**

Vacuum Pulse Setpoint **20.0** inHG

Drain Phase Time **1** Min.

Drying Phase Time **2** Min.

CV Enabled in Drying **YES**

Cycle Start

- Previous Picture
- Splash Screen
- Synoptic
- Freeze Dry Screen
- Leak Test Screen
- Function Test Screen
- Sterilization
- CIP
- CIP Skid
- Defrost
- Semi-Auto Functions
- Historical Trend
- Product Probe Screen
- Maintenance Screen
- Cycle Continue
- End Freeze Dry Cycle
- Cancel Current Cycle
- Historical Export
- Print Screen

Cycle Status


Cycle	Clean in Place
Total Time	0 Min.
Phase	Ready
Phase Time	0 Min.
Step	0
Step Time	0 Min.
Vacuum	1025 mT
Shelf Temp.	20.3 °C
Cond. Temp.	17.7 °C

Alarm Summary



SUPER


CIP Skid



Benchmark Control Platforms


Encore

Allen Bradley Micrologix controller. Many of the features available on the more expensive Maestro but limited options to keep the cost down. BE, EC-1, CM, PRCM security historical data export all available (Standard)



Maestro

SLC 5/03 AB controls offers world wide support for parts and service. The Maestro has additional options such as 21cfr11 compliant security and System Integration Testing, unlimited product probes, customized programming. Best choice for high tech R& D facilities or GMP production.




Controls

Basic Components:


White – Vacuum

Yellow – Heat Transfer



Blue – Refrigeration



Redundant Compressors and pumps



Full N + 1 Redundancy



Intuitive Recipe Screen

Freeze Dry Recipe

Cycle Name **SPMASTER**
 Batch ID **Herold**
 Shelf Load Temperature **26.0** Deg C

05/25/2004 08:07:48

Shelf Load

Thermal Treatment

Step	Rate/Hold	Temperature	Time
1	hold	0.0 Deg C	0 Min.
2	hold	0.0 Deg C	0 Min.
3	hold	0.0 Deg C	0 Min.
4	hold	0.0 Deg C	0 Min.
5	hold	0.0 Deg C	0 Min.
6	hold	0.0 Deg C	0 Min.
7	hold	0.0 Deg C	0 Min.
8	hold	0.0 Deg C	0 Min.
9	hold	0.0 Deg C	0 Min.
10	hold	0.0 Deg C	0 Min.
11	hold	0.0 Deg C	0 Min.
12	hold	0.0 Deg C	0 Min.

Freeze, Condenser and Evacuate

Freeze Temperature **25.0** Deg C
 Extra Freeze Time **0** Min.
 Vacuum Start Permit (Condenser Temp.) **-20.0** Deg C
 Heat Start Permit (Vacuum) **700** mT

Primary Drying

Step	Rate/Hold	Temperature	Time	Pressure	Press. Rise
1	hold	25.0 Deg C	1 Min.	600 mT	0 mT
2	hold	0.0 Deg C	0 Min.	0 mT	0 mT
3	hold	0.0 Deg C	0 Min.	0 mT	0 mT
4	hold	0.0 Deg C	0 Min.	0 mT	0 mT
5	hold	0.0 Deg C	0 Min.	0 mT	0 mT
6	hold	0.0 Deg C	0 Min.	0 mT	0 mT
7	hold	0.0 Deg C	0 Min.	0 mT	0 mT
8	hold	0.0 Deg C	0 Min.	0 mT	0 mT
9	hold	0.0 Deg C	0 Min.	0 mT	0 mT
10	hold	0.0 Deg C	0 Min.	0 mT	0 mT
11	hold	0.0 Deg C	0 Min.	0 mT	0 mT
12	hold	0.0 Deg C	0 Min.	0 mT	0 mT
13	hold	0.0 Deg C	0 Min.	0 mT	0 mT
14	hold	0.0 Deg C	0 Min.	0 mT	0 mT
15	hold	0.0 Deg C	0 Min.	0 mT	0 mT
16	hold	0.0 Deg C	0 Min.	0 mT	0 mT

Pressure Rise Test **NO** Pressure Control Action **no**
 Valve Closed Time **5** Sec. Repeat Test Time **5** Min.

Primary Instant P_Rlse

Secondary Drying

Step	Rate/Hold	Temperature	Time	Pressure	Press. Rise
1	hold	24.0 Deg C	1 Min.	550 mT	0 mT

Product Temperature **0.0** Deg C

Pressure Rise Test **NO** Pressure Control Action **no**
 Valve Closed Time **5** Sec. Repeat Test Time **5** Min.
 Delay Start Test Time **5** Min.

Secondary Instant P_Rlse

Storage Temperature **23.0** Deg C

Previous Picture

Splash Screen

Synoptic

Freeze Dry Screen

Leak Test Screen

Function Test Screen

Sterilization

CIP

Defrost

Semi-Auto Functions

Historical Trend

Product Probe Screen

Maintenance Screen

Cycle Continue

End Freeze Dry Cycle

Cancel Current Cycle

Historical Export

Print Screen

Cycle Status

Cycle Freeze Drying
 Total Time 0 Min.
 Phase Ready
 Phase Time 0 Min.
 Step 0
 Step Time 0 Min.
 Vacuum 10000 mT
 Shelf Temp. 24.0 °C
 Cond. Temp. 23.4 °C

Download Recipe

Recipe Status

Completed - Okay

Start

Recipe Manager

Recipe Modification


Alarm Summary

SUPERVISOR



Powerful Historical Trend Screen

Allows data to be viewed in any interval for any available recorded data.



Manual Operation for Easy Troubleshooting

05/07/2004 13:26:23

Semi-Auto

Freeze Condenser Vacuum Heat Vacuum Release Isolation Valve Drain Water Ring Pump Component Mode

250 mT 45.0 °C

- Previous Picture
- Splash Screen
- Synoptic
- Freeze Dry Screen
- Leak Test Screen
- Function Test Screen
- Sterilization
- CIP
- Defrost
- Semi-Auto Functions
- Historical Trend
- Product Probe Screen
- Maintenance Screen
- Cycle Continue
- End Freeze Dry Cycle
- Cancel Current Cycle
- Historical Export
- Print Screen

Cycle Status

Cycle	Steam Defrost
Total Time	0 Min.
Phase	Ready
Phase Time	0 Min.
Step	0
Step Time	0 Min.
Vacuum	10000 mT
Shelf Temp.	24.4 °C
Cond. Temp.	-7.6 °C

Alarm Summary SUPERVISOR



Controls Considerations

Recipe requirements?

How much time do you have to get it right?

Process development tools:

Smart
BE or PVG/CM (See Tech Briefs)

Automated Testing

Secure Data

21CFR11?

Freeze Dry Recipe

05/25/2004 08:07:48
 Cycle Name: SPMMASTER
 Batch ID: Herold
 Shelf Load Temperature: 26.0 Deg C

Thermal Treatment

Step	Rate/Hold	Temperature	Time
1	hold	0.0 Deg C	0 Min.
2	hold	0.0 Deg C	0 Min.
3	hold	0.0 Deg C	0 Min.
4	hold	0.0 Deg C	0 Min.
5	hold	0.0 Deg C	0 Min.
6	hold	0.0 Deg C	0 Min.
7	hold	0.0 Deg C	0 Min.
8	hold	0.0 Deg C	0 Min.
9	hold	0.0 Deg C	0 Min.
10	hold	0.0 Deg C	0 Min.
11	hold	0.0 Deg C	0 Min.
12	hold	0.0 Deg C	0 Min.

Primary Drying

Step	Rate/Hold	Temperature	Time	Pressure	Press. Rise
1	hold	25.0 Deg C	0 Min.	600 mT	0 mT
2	hold	0.0 Deg C	0 Min.	0 mT	0 mT
3	hold	0.0 Deg C	0 Min.	0 mT	0 mT
4	hold	0.0 Deg C	0 Min.	0 mT	0 mT
5	hold	0.0 Deg C	0 Min.	0 mT	0 mT
6	hold	0.0 Deg C	0 Min.	0 mT	0 mT
7	hold	0.0 Deg C	0 Min.	0 mT	0 mT
8	hold	0.0 Deg C	0 Min.	0 mT	0 mT
9	hold	0.0 Deg C	0 Min.	0 mT	0 mT
10	hold	0.0 Deg C	0 Min.	0 mT	0 mT
11	hold	0.0 Deg C	0 Min.	0 mT	0 mT
12	hold	0.0 Deg C	0 Min.	0 mT	0 mT
13	hold	0.0 Deg C	0 Min.	0 mT	0 mT
14	hold	0.0 Deg C	0 Min.	0 mT	0 mT
15	hold	0.0 Deg C	0 Min.	0 mT	0 mT
16	hold	0.0 Deg C	0 Min.	0 mT	0 mT

Freeze, Condenser and Evacuate

Freeze Temperature: 25.0 Deg C
 Extra Freeze Time: 0 Min.
 Vacuum Start Permit (Condenser Temp): 20.0 Deg C
 Heat Start Permit (Vacuum): 700 mT

Pressure Rise Test: NO Pressure Control Action: no
 Valve Closed Time: 5 Sec. Repeat Test Time: 5 Min.

Secondary Drying

Step	Rate/Hold	Temperature	Time	Pressure	Press. Rise
1	hold	24.0 Deg C	1 Min.	550 mT	0 mT

Product Temperature: 0.0 Deg C
 Pressure Rise Test: NO Pressure Control Action: no
 Valve Closed Time: 5 Sec. Repeat Test Time: 5 Min.
 Delay Start Test Time: 5 Min.

Storage Temperature: 23.0 Deg C

Download Recipe
 Recipe Status: Completed - Okay

Start Recipe Manager Recipe Modification

Alarm Summary SUPERVISOR

Cycle Status

Cycle: Freeze Drying
 Total Time: 0 Min.
 Phase: Ready
 Phase Time: 0 Min.
 Step: 0
 Step Time: 0 Min.
 Vacuum: 10000 mT
 Shelf Temp: 24.0 °C
 Cond. Temp: 23.4 °C




Recipe Screen




Historical Data

Pirani (T/C) Vacuum Gauge

- Wire heated by electric current
- Molecules of process gas impinge upon wire and cool it
- More molecules = cooler
- Less molecules = hotter
- Temperature correlates to Pressure
- Must be calibrated for specific process gas




Capacitance Manometer



Encore


Cycle Name: CycleName
Batch ID: batchid

Temperature/Vacuum



Condenser Temperature # 1	-84.4
Shelf Temperature	10.1
Temperature Reference (shelf)	10.0
Vacuum [CM]	100.0
Vacuum [PVG]	103.0
Vacuum Reference	100.0

Agreement with PVG/CM and Probes



Electrical Cabinet



Skid Mounted UPS for PLC & Workstation



Documentation & Validation

Factory and Site Acceptance Testing (FAT/SAT)


Customized Workbooks with Material and Welding Certifications

Control System Integration Testing

Inverted GAMP 4 Triangle

User Requirement Specs
Functional and Design Specs

Operational Qualification



Service & Support

Installation Supervision

Start-up Supervision

**Maintenance, Operator or
Lyophilization Training**

PM Contracts

Extended Warrantee

**Emergency Weekend
Response**

Time and Materials

Maintenance Audit Programs

**Calibration / Validation
Assistance**





VirTis



SP SCIENTIFIC

Partnership Through Support ...

... Together we can get the
job done better!



Service and Support Mission

Our mission is Clear ... to provide support and services of the highest quality, customer experience and greatest possible value to our customers, thereby gaining an holding their respect and loyalty.

In essence Partnership through Support

