





# FEATURES

- Low service costs
- Long life
- Perfect for storage of rubber and steel for testing



# PRO

## LOW TEMPERATURE FOR INDUSTRIAL RESEARCH



PRO 11



PRO 21



PRO 31



PRO 41



PRO 60

	PRO 11	PRO 21	PRO 31	PRO 41	PRO 60
Ambient temperature (Celsius)	+10 to +30	+10 to +30	+10 to +30	+10 to +30	+10 to +30
Gross volume, litres	130	226	300	360	416
Net volume, litres	130	226	300	360	416
Temperature range inside unit	-30 to -60	-30 to -60	-30 to -60	-30 to -60	-30 to -60
Refrigerant: Mix of 2 Hydro carbons	Yes	Yes	Yes	Yes	Yes
Input power (watt):	375	455	544	587	612
Energy consumption / 24 h (kwh) 50 % running time	4,3	5	5,2	6,1	9,2
<b>DIMENSIONS</b>					
External dimensions (wxdxh)	725*655*865	1055*655*865	1300*655*865	1500*655*865	1700*655*865
Internal dimensions (wxdxh)	520*450*650	850*450*650	1100*450*650	1300*450*650	1500*450*650
Dimension A	725	1055	1300	1500	1700
Dimension B	655	655	655	655	655
Dimension C	865	865	865	865	865
Dimension D	780	780	780	780	780
Dimension E	520	850	1100	1300	1500
Dimension F	320	650	900	1100	1300
Dimension G	450	450	450	450	450
Dimension H	650	650	650	650	650
Dimension J	400	400	400	400	400
Dimension L	850	850	850	850	850
Insulation ,mm	100	100	100	100	80
<b>ACCESSORIES</b>					
Dividers	Optional	Optional	Optional	Optional	Optional
Max. number of baskets (option)	1	2	3	4	4
Castors	Optional	Optional	Optional	Optional	Optional
Castors w. brakes	Optional	Optional	Optional	Optional	Optional
Lock	Included	Included	Included	Included	Included
Potential free contact included	Yes	Yes	Yes	Yes	Yes
Digital controller w. acoustic alarm	Included	Included	Included	Included	Included
<b>OTHER DATA</b>					
Outer cabinet	White	White	White	White	White
Inner cabinet	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium
Type of cooling	Static	Static	Static	Static	Static
Type of defrost	Manual	Manual	Manual	Manual	Manual
Power	220V/50Hz	220V/50hz	220V/50hz	220V/50hz	220V/50hz
Power 110volt / 60 hz	Optional	Optional	Optional	Optional	Optional
Thermostat	Digital	Digital	Digital	Digital	Digital
Net weight, kg	56	67	80	89	95
Gross weight, kg	63	74	88	97	104
<b>SHIPPING</b>					
20' cont. Capacity	42	30	24	20	18
40' cont. Capacity	90	64	50	44	38
Standard Trailer capacity	102	72	56	50	44