



Autoclave for sterilization “MED 12”

DESKTOP MODEL

Features: see page 58



MODEL

Part No.	Capacity litres	Ø (usable) cm	Height (usable) cm	Height/Width/Depth (exterior) cm	Power W	Weight Kg
4001756	12	24	26	46 38 38	1800	18

It comes complete with a stainless steel wire basket of Ø 22x24cm high. Part no. 1004772.



Autoclave for sterilization “MED 20”

VERTICAL STAND MODEL.

CONFORMS TO APPARATUS PRESSURE REGULATION DIRECTIVES.

Features: see page 58



MODEL

Part No.	Capacity litres	Ø (usable) cm	Height (usable) cm	Height/Width/Depth (exterior) cm	Power W	Weight Kg
4001757	20	24	41	70 38 38	2500	41

It comes complete with two stainless steel wire baskets of Ø 22x24cm and Ø 22x14cm high. Part no. 1004772 and 1004771.