



Minifors

Unpack, Install and Ferment



- ▶ Designed for easy handling
- ▶ Simple Operation
- ▶ Complete Packages
- ▶ Excellent Value
- ▶ Builds Bridges



www.infors-ht.com

May be subject to technical amendments

INFORS HT

Maximum capability for minimum price

▶ Designed for Easy Handling

Minifors is designed to make handling as easy as possible for new users with „plug'n'pump“, the use of a heating block rather than a jacket and an integral frame for carrying all components to the autoclave together.

▶ Simple Operation

Setting up the Minifors is made as easy as possible with a control panel that is very similar to that used on our Multitron incubator shakers.

▶ Complete Packages

Minifors comes as a package with Cell and Bacterial versions plus two vessel sizes. It can be set up „out of the box“ with the necessary accessories provided.

▶ Excellent Value

The price/performance ratio of Minifors is high with a full system being comparable in cost to a large incubator shaker but with the ability to provide more culture in less space.

▶ Builds Bridges

Minifors provides the link between shake flasks and larger-scale fermentation by providing the shake flask user with the additional features they need most i.e. volume, sampling, pH measurement/control, oxygen control, the ability to add feed and PC data-logging/control.



Complete Specifications Additional parameters

The basic system delivers practically all the usual process parameters for fermentation.

- Stirrer speed
- Temperature
- pH
- pO₂
- Antifoam/level
- Feed
- PC Software: Iris

With the optional 2 additional analogue, inputs/outputs, the system can cater for additional parameters if required, such as, Mass-Flow, turbidity, exit gas analysis and additional external analogue pumps.

Infors AG
Headoffice, Switzerland

Rittergasse 27
CH-4103 Bottmingen
T +41 (0)61 425 77 00
F +41 (0)61 425 77 01
headoffice@infors-ht.com

Infors GmbH

Dachauer Str. 6
D-85254 Einsbach
T +49 (0)8135 8333
F +49 (0)8135 8320
infors.de@infors-ht.com

Infors UK Ltd

The Courtyard Business
Centre Dovers Farm,
Lonesome Lane Reigate,
Surrey, RH2 7QT, UK
T +44 (0)1737 22 31 00
F +44 (0)1737 24 72 13
infors.uk@infors-ht.com

Infors Sarl

2 rue du Buisson
aux Fraises
Bâtiment D13
FR-91300 Massy
T +33 (0)1 69 30 95 04
F +33 (0)1 69 30 95 05
infors.fr@infors-ht.com

Infors Benelux BV

Fabriekstraat 38
NL-7005 AR Doetinchem
P.O. Box 613
NL-7000 AP Doetinchem
T +31 (0)314 364 450
F +31 (0)314 378 076
infors.bnl@infors-ht.com

Infors Canada

8350 rue Bombardier
Anjou, Quebec
Canada H1J 1A6
T + 1 514 352 50 95
F + 1 514 352 56 10
infors.ca@infors-ht.com

Infors Biotechnology
(China) Co., Ltd.
Room 505C, Building 106
Lize Zhongyuan
Wangjing New Industrial Zone
Chaoyang District, Beijing
100102 P. R. of China
T +86 (0)10 51 65 20 68
F +86 (0)10 64 39 05 85
info@infors-ht.com.cn