# VirTis Freezemobile 35L

Research Freeze Dryer



### **Electrical Requirements**

#### XL EL Item Voltage (VAC) 208 / 240 208 / 240 50, 60 50, 60 Hertz Phase 1 1 Breaker Amperage 30 20 **NEMA Receptacle** 6-20R\* L6-30R\*

### **Utility Requirements**

Item	XL	EL
Approx. Peak Heat Generated (BTU/h)	7,200	9,800

#### **Key Features**

- Compact, freestanding, mobile design.
- Includes microprocessor-based Sentry 2.0 controller.
- Choice of refrigeration system to meet a variety of process requirements.
- Wide variety of optional manifolds, shelves and accessories available.
- Optional filter traps available to protect the vacuum pump from corrosive solvents.

### **Performance Specifications**

Item	XL	EL
Lowest Condenser Temperature (50 Hz / 60 Hz, °C)	-67 / -70	-82 / -85
Lowest Freezing Temperature with Optional Shell Bath (50 Hz / 60 Hz, °C)	-52 / -55	-67 / -70
Maximum Condenser Capacity (L)	35	35
Condenser Surface Area (in <sup>2</sup> / cm <sup>2</sup> )	750 / 4838	750 / 4838
Maximum Ice Condensing Capacity in 24 hours (L) <sup>†</sup>	20	20
Maximum Deposition Rate (L/hour) <sup>†</sup>	0.83	0.83
Condenser Pull-Down from 20 °C to -45 °C (minutes)	≤ 25	≤ 25
Number of Compressors	1	2
Compressor Horsepower	1.5	1.5, 1.0
System Refrigerant	R245fa / R508B	R508B, R407C
Vacuum Time to 100 Millitorr (minutes)	≤ 20	≤ 20
Lowest System Vacuum (mT)	≤ 15	≤ 15
Average Defrost Time (minutes)	≤ 60	≤ 60

**Note**: Performance specifications are based on SP Scientific test data from units operating at an ambient room temperature of approximately 20 °C. SP Industries recommends an operating range of 15-25 °C (59-77 °F).



## VirTis Freezemobile 35L

### Research Freeze Dryer





### **Dimensional Data**

Width	35 in (89 cm)
Depth	36 in (91 cm)
Height	37 in (94 cm)
Approximate Weight Range	450-500 lb (205-225 kg)

### **Additional Information**

Construction	316L Stainless Steel Condenser Chamber
Vacuum Pump (required, not included)	Two-Stage Rotary Vane
Defrost Type	Hot Gas
Refrigerant	CFC-Free

Assembly

4, 8, 12, 24

or 48 Port

Stainless

Steel Vertical

Manifold



Acrylic Drum Manifold

12 Port Acrylic Drum Manifold; 12 Inch Outside Diameter (OD)



**Stainless Steel Drum** Manifold

12 or 18 Port Stainless Steel Drum Manifold

### **Optional Components**



**Bulk Shelf** Rack

3, 4 or 5 Shelves: Heated or Unheated



Stoppering Vertical Manifold Assembly

1, 2 or 3 Shelves; Heated or Unheated



Horizontal "T" Manifold

> 1, 2 or 3 Tier Stainless Steel Horizontal "T' Manifold

Note: Additional standard and custom accessories are available, along with flask adapters, glassware and other components. Contact SP Scientific for more information.

\* Receptacle configuration for 50 Hz units varies by country.

<sup>†</sup> The specified Maximum Ice Condensing Capacity in 24 Hours and Maximum Deposition Rate are based on the process of freeze-drying water as aggressively as possible. The freeze dryer's ability to collect ice at an hourly rate or over a specified period will always be application dependent.

Copyright © 2013 SP Scientific. All marks herein are used under license. All brand or product names mentioned may be trademarks or registered trademarks of their respective companies. SP Scientific reserves the right to change specifications without notice.