

T-CLASSIC

9-Litre Portable Autoclave

Making your process easy, quick and affordable.

The T-Classic One, is a lightweight, portable steam autoclave designed as a cost-effective sterilization solution. Fully automatic with an easy one touch operation, delivering short cycles.

Ideally designed to meet portability needs of a beauty clinic or academic facility.

The T-Classic One can also be used in combination with another autoclave or as a back up.

Tuttnauer Performance with the simplicity and affordability your small practice needs



T-Classic One

Instruments sterilized and ready quickly



The T-Classic One offers competitively quick unwrapped cycle times, in as short as 18 minutes*. The fast point-of-use sterilization makes your turnaround a priority, increasing efficiency, while eliminating bottlenecks.

*Cycle times may vary, depending on load & weight.



The T-Classic One is a simple one step operated device. With one touch, the T-Classic One will begin an automatic cycle until completion. An electronic detector will turn off the power if water level is insufficient. The perfect autoclave for limited space.



The T-Classic One is the lightest tabletop autoclave available, weighing only 5.35 kg allowing for easy portability, while requiring minimal counter space.